

#15

Express Mail Label No.: EL326923155US
Date of Deposit: November 9, 2001

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Application of:	Carstens, Carsten-Peter	Examiner:	G. Leffers Jr.
Serial No.:	09/492,590	Group:	1636
Filed:	January 27, 2000		
Entitled:	METHODS AND COMPOSITIONS FOR HIGH LEVEL EXPRESSION OF A HETEROLOGOUS PROTEIN WITH POOR CODON USAGE		

Attorney Docket No.: 25436/1340 (Formerly 4114/85530)

Commissioner for Patents
Washington, D.C. 20231

DECLARATION UNDER 37 C.F.R. §1.131

I declare:

1. I, Carsten-Peter Carstens, am the inventor of the invention claimed in the above-noted U.S. patent application.
2. I have read and understood the Office Action mailed May 9, 2001, and have read and understood the cited reference, U.S. Patent No. 6,214,602 ("the '602 patent," issued to Zdanovsky et al. on April 10, 2001, from an application filed August 28, 1998). I understand that the Examiner has cited the '602 patent as a novelty reference over claims 1-10, 15, 16, 19, 22-27, 32-40 and 42-44, and as an obviousness reference over claims 18, 20 and 21. The '602 patent is cited as teaching vectors and methods for the overproduction of Clostridium toxins and proteins by hosts such as E. coli. The methods taught by the '602 patent are said to feature the use of host cells containing a recombinant expression vector, wherein the expression vector encodes tRNAs that recognize rare codons and wherein the host cell expresses at least a fragment of at least one clostridial protein. The patent is also said to teach vector constructs that encode three tRNAs that recognize rarely used codons.
3. Prior to the August 28, 1998 filing date of the '602 patent, I had both conceived of and reduced to practice the claimed invention. The attached exhibits A and B, consisting of copies of my notebook entries detail the experiments that gave rise to the claimed invention, all performed before August 28, 1998. The dates on these exhibits have been redacted. For convenience, Appendices I and II contain transcripts of the relevant notebook entries.

A. Vector comprising an array of three or more tRNA genes which correspond to codons that are rarely used in a host cell

The first set of notebook entries, Exhibit A (transcripts in Appendix I), describes the development of pACYC-based vectors encoding genes for three recombinant tRNAs specific for rarely used codons. The first array is termed "RIL," for *ArgU* (R), *IleY* (I) and *LeuW* (L).

First, even the initial entry of Exhibit A (Entry # 1), describing experiments performed before August 28, 1998, shows that the vectors comprising an array of three tRNA genes were conceived of even before that date. That is, because Entry # 1 details efforts to ligate the RIL insert into the pACYC-LIC vector, I already had the critical RIL insert fragment and the intent to put it into an expression vector that would meet the limitations of the claims.

Second, Entry # 6 shows the ligation that generated pACYC-RIL vector constructs 49, 50 and 55, which were shown in the ensuing experiments (e.g., the PCR screen of Entry # 8 and the miniprep digests of Entry # 9) to contain the RIL insert, all performed before August 28, 1998. Thus, the vectors comprising an array of three or more tRNA genes which correspond to rarely used codons as claimed in the present application were reduced to practice as of the date of entry of Entry # 6, which is prior to August 28, 1998. The vectors resulting from the ligation were screened, isolated and confirmed by the date of Entry # 9, which is also before August 28, 1998. While it is clear that the vectors were made as of the date of Entry # 6 but not confirmed until the date of Entry # 9, both experiments were performed in advance of the August 28, 1998 priority date of the '602 patent.

B. Host cell containing a recombinant DNA molecule which comprises an array of three or more tRNA genes, wherein the tRNA genes correspond to codons that are rarely used in the host cell

The intent at the time of making the pACYC-RIL constructs was to introduce them to host cells for their influence on the expression of heterologous genes with rare codons. Thus, the conception of the claimed host cells that comprise an array of three or more tRNA genes, wherein the tRNA genes correspond to codons that are rarely used in the host cell, necessarily preceded, or at a minimum, coincided with, the conception of the vectors themselves. Therefore, Exhibit A and Appendix I are also sufficient to demonstrate conception of the claimed host cells before August 28, 1998.

Similarly, Exhibit A and Appendix I are sufficient to demonstrate reduction to practice of host cells comprising an array of three or more tRNA genes, wherein the tRNA genes correspond to codons that are rarely used in the host cells. Specifically, the transformation of XL1-Blue host cells with the ligation reactions of Entry # 6 created host cells carrying such a construct. The PCR screening documented in Entry # 7 and the miniprep digests of Entry # 9 show that the pACYC-RIL clones 49, 50 and 55 are positive for the RIL insert. Finally, the test of Fip expression documented in Entries 10, 11 and 12 shows the introduction and induction of the pACYC-RIL construct in the well-

known protein expression *E. coli* host strain BL21(DE3). Each of these experiments was carried out before August 28, 1998.

In view of the above, host cells as claimed were conceived of and reduced to practice in advance of the August 28, 1998 priority date of the '002 patent.

C. Method of producing a protein of interest

The second set of notebook entries, Exhibit B (transcripts in Appendix II), describes the use of the pACYC-RIL constructs to increase the production of a protein of interest from a gene with a non-*E. coli* codon bias.

The conception of a method of producing a protein from a gene with infrequently used codons, the method involving the expression of specific recombinant tRNAs from a vector in a host cell is necessarily very closely linked in time, if not coincident with or even prior to, the conception of the vector and host. Thus, conception of the method occurred before the performing the experiments described in the earliest entry in Exhibit A, which is prior to August 28, 1998. The reduction to practice is described below.

The hyperthermophilic archaeon *Pyrococcus furiosus* has a very AT-rich genome that results in distinctly different codon usage from the less AT-rich *E. coli* host cell genome. Therefore, one wishing to express *Pfu* DNA polymerase in *E. coli* will achieve only limited expression unless codon bias is compensated.

Exhibit B notebook Entry # 1 (page 48) shows the transformation of BL21(DE3)-RIL cells (corresponds to isolate 49) with a *Pfu* polymerase fusion protein. Entry # 2 shows a Coomassie-stained gel containing host cell proteins from *Pfu* polymerase construct-transformed cells with and without IPTG induction of tRNA construct expression. The gel photo shows the production of the *Pfu* polymerase fusion protein. The gel also shows, and the notebook entry states, that the production of the *Pfu* polymerase (rich in rare arg codons and Ile codons), is enhanced in RIL 49 host cells. Both experiments have been performed before August 28, 1998.

Exhibit B Entry # 2 therefore shows reduction to practice for a method of producing a protein of interest, comprising the step of culturing a host cell containing a recombinant DNA molecule that comprises an array of three or more tRNA genes, wherein the tRNA genes correspond to codons that are rarely used in the host cell, wherein the codons are present in the gene for the protein of interest, and wherein the conditions of culturing the host cell are sufficient to produce the protein of interest. The reduction to practice thus occurred prior to August 28, 1998.

4. I hereby declare that all statements made herein of my own knowledge are true, and that all statements made on information and belief are believed to be true; and further, that these statements were made with the knowledge that willful, false statements and the like so made are punishable by fine or imprisonment, or both, under Section 1001 of Title 18 of the United States Code, and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

11/18/01
Date


Carsten-Peter Carstens

Appendix I

Transcripts of notebook entries relating to vectors and host cells

(Note -- several unrelated constructs were also being made in addition to the tRNA array constructs):

Entry #1

[page 066]

annealing of LIC-cloning reactions

1 μ l of LIC treated vector (pACYC-LIC or pESP4)
5 μ l LIC-treated inserts (Flp, Ku70 and Ku80 and RIL)
4 μ l dH₂O
incubated for 1 h at RT

vectors setup: pESP4Ku70, pESP4Ku80, pESP4Flp and pACYC-Flp and RIL

used 1 μ l to transform XL1-Blue cells. Plated transformations on LB-Amp (pESP4 clones) or LB-Cam plates (pACYC based).

Entry #2

[page 067]

PCR Screening of pACYCFlp, pACYC-RIL, pESP4-Ku70, pESP4-Ku80 and pESP4-Flp colonies (Continuation of previous page)

[details PCR reaction components and cycling conditions]

[bottom of page 067, below gel photo]

None of the pACYC-RIL screened are positive.

Entry #3

[page 068]

Screening of more pACYC RIL clones (cont' from prev. page)

[details reaction setup and cycling conditions and shows diagnostic gel; no positive clones detected]

Confirmation digests of RIL fragment used for LIC cloning of pACYC-RIL and miniprep digest of selected ESP4Flp, ESP4 Ku80 and ESP4 Ku70 clones

[bottom of page, next to lower gel photo]

test of RIL fragment

undigested: 610 bp

NcoI: 220 bp, 390 bp (high salt retards small frag. -- NcoI uses 2X UB
(universal buffer))

XbaI: 180, 430 bp (+LIC sites)

the RIL fragment seems to be correct
unclear why I can't clone - toxic?

Entry #4

[page 080]

Ordered new primer for construction of RIL array. LIC cloning appeared to be very ineffective in previous attempts. Try to clone between the SpeI and XhoI sites of pACYC-LIC. Also ordered forward and reverse primer for RI and IL array to be cloned the same way.

Entry #5

[page 084]

PCR-amplification of RIL, IL and RI fragments using primer with Spe/Xho extensions

[details reaction setup and cycling conditions and isolation of amplified fragments following gel electrophoresis]

Entry #6

[page 086]

Construction of pACYC-RIL, pACYC-RI and pACYC-IL.

- a) vector digest [details digest of pACYC-LIC with SpeI and XhoI]
- b) ligation [details ligation using 2 μ l of cut pACYC-LIC vector and 5 ml of RIL, RI or IL]

used 2 μ l to transform XL1 Blue

Entry #7

[page 087]

PCR screen of pACYC-RIL and pACYC-IL (see prev. page)

[details reaction setup and cycling conditions and shows gel photo]

expected sizes:

RIL 610 bp

IL 430 bp

pACYC-RIL: clone #7

pACYC-RI: 23, 25, 28, 29, 32, 33

Entry #8

[page 090]

minipreps of pACYC-RI and pACYC-RIL colonies (continuation of page 87)

expected fragment sizes

2192, 1230 (RI) or 1400 (RIL)

both pACYC-RI clones are confirmed by restriction digests.

None of the RIL clones is positive although they were PCR positive (comp. page 87)

[page 091]

Rescreened pACYC RIL colonies by PCR

[details reaction setup and cycling conditions, shows gel photo]

none of the previous picked colonies are positive.

3 of the new colonies (49, 50 and 55) are positive.

Entry #9

[page 091]

minipreps of pACYC-RIL clones (above)

expected

NcoI: 2192 bp, 1400 bp

XbaI:

All clones are positive by NcoI digests, but XbaI site seems to be missing (either from vector or from insert).

Entry #10

[page 095]

test of Flp expression in pACYC-RIL/BL21DE3 cells

picked colonies from transformations on page 93 (one each expression clone in each RIL strain) and inoculated O/N cultures (1 ml Amp/Cam LB)

Entry #11

[page 095]

induction of Flp expression in BL21DE3/RIL strains

[details growth and induction of expression of Flp]

Entry #12

[page 096]

detection of Flp expression in BL21DE3 RIL strains

[details gradient gel electrophoresis and Western blotting of samples from cells carrying RIL 49, RIL 50 and RIL 55 constructs]

Signal extremely strong. Had to add piece of paper between the membrane and the film to achieve reasonable exposure times (chemilum inogram next page). Took image with Eagle eye (signal strong enough for that).

Appendix II

Transcripts of notebook entries relating to production of a protein of interest

Entry #1

[page 48]

Comp IRL cells

[details making BL21DE3-IRL cells competent for transformation]

transformed 100 µl ea. with

Pfu polymerase (CBP fusion)	} archaeal protein
Flavo PEF	} AT rich
CreProL#5	

plated on LB Amp.

Entry #2

[page 55]

Coomassie gel of NusA production and IRL test.

production of Pfu-polymerase

(rich in rare arg codons and IleY codons) is enhanced in IRL49. At least the argU component of the strain enhances Pfu production.



inoculated cells ~~from~~ ^{Be An} for #5, 6 and 4 for miniamp

inoculated on culture of BL21 DEBIL for comp cells

comp RIL cells

inoculated 10 ml LB with 100 µl O/L culture

↓
let grow for 3h at 37°C in shaker

↓
pellet cells for 15 at 2500 rpm in table top

↓
resuspension 2.5 ml F&B

↓
pellet 15 at 2500 rpm

↓
resuspended in 600 µl F&B

↓
transformed 100 µl cell. with

Plu Polymerase (C&G fusion) } extract protein
Havus PEF } At risk
Cre Prot #5

plated on LB Amp

Confirmation PCR of clones 56 well (miniprep DNA)

	A	B	C	
10xTaq poly	30	30	30	95°C, 3' / 55°C, 3' / 72°C, 2' / 10 cycles
10xTaq poly b.	30	30	30	95°C, 1' / 55°C, 1' / 72°C, 2' / 24 cycles
AMPs (25 min. eq)	2.4	2.4	2.4	72°C, 10' → 60°C, hold
1047	.5	-	-	
1048	-	.5	-	loaded 1 µl each sample on
51C	.5	.5	-	10% agarose gel
53C	-	-	.5	
54C	-	-	.5	image next page
PEF	-	-	-	
Taq 2020	2	2	-	

Q. Q. Q.



EXHIBIT A

066

ENTRY # 1

annealing of LIC-cloning reactions

- 1 μ l LIC-treated vector (pACYC-LIC or pESP4)
- 5 μ l LIC-treated inserts (Fbp, ku20 and ku80 and RIL)
- 4 μ l ddH₂O
- incubated for 1h at RT

Gelatin setup: pESP4 ku20, pESP4 ku80, pESP4 Fbp and pACYC-Fbp and RIL

goal of pACYC-Fbp \rightarrow Fbp activity is easier to detect in an in vitro assay \rightarrow quality control of the insert

used 1 μ l to transform XL-blue cells. Plated transformations on LB Amp (ESP4 clones) or LB Can plates (pACYC based)

PCR amplification of kanamycin resistance marker from pACYC-177
for confirmation of correct deletion of pESP4

	1	2
10x Pfu pol buffer (cloned)	10	20
dNTPs (25mM)	.8	.6
PCR	1	1
pACYC-177	1	1
CC 1830 (kan fwd +202)	.5	-
CC 1834 (kan rev)	-	.5
CC 1156 (kan rev)	.5	.5
cloned Pfu pol	2	2
ddH ₂ O	84.2	84.2

used hot top

95°C, 3' / 55°C, 3' / 72°C, 3' // 1 cycle


95°C, 1' / 55°C, 1' / 72°C, 3' // 24 cycles

72°C, 10' \rightarrow 6°C, ∞

loaded reactions on 1% agarose gel

NE EAGLE EYE II 04/03/98 17:36:46

26 COUNTS.
17:36:16 1998.



Expected product size: 1.1 kb isolated fragments using Strataprep gel extraction kit
eluted DNA into 50 μ l TE


[Signature]

PCR screening of pACYC-Flp, pACYC-RIL, pESP4-Ku70, pESP4-Ku80 and pESP4-Flp
colonies (Continuation of previous page)

	1	2	3	4	
10x Tuff poly buffer	25	25	75	100	used 47.5 μ l / reaction
dNTPs (25mM)	2	2	6	8	
CC 937 (Ku70 screen)	.5	-	-	-	95°C, 3' / 55°C, 3' / 72°C, 2' // 1 cycle
CC 938 (Ku80 screen)	.5	-	-	-	95°C, 2' / 55°C, 1' / 72°C, 2' // 29 cycles
CC 939 } Ku80 screen	-	.5	-	-	72°C, 10' \rightarrow hold at 6°C
CC 940 } Ku80 screen	-	.5	-	-	used hot nap
CC 1375 } Flp	-	-	.5	-	
CC 1376 } Flp	-	-	.5	-	loaded 10 μ l each reaction on
CC 1605 (control band)	-	-	-	1	1% agarose gel
CC 1725 (Ku70 rev.)	-	-	-	1	
dist Tuff 2000	25	25	75	10	
dist H ₂ O	214.5	214.5	660	880	

STRATAGENE EAGLE EYE II 04/06/98 16:09:59

IMAGE SIZE (640 x 480 x 8).
INT PERIOD = 0.16 SEC.
ACQUIRED MON APR 06 16:08:46 1998.



expected fragment product sizes

Ku70:

Ku80:

Flp : 1270 bp

RIL : 610 bp

all Ku70 clones, the
Ku80 clones 6 and 7,
ESP4 Flp clones 11 and 13-15
and pACYC clone 23 are
positive in PCR screen
 \rightarrow product for subcloning DNA


none of the pACYC-RIL clones
are positive

Xing-Hu

C.C.D.

Screening of more pACYC-RIL clones (cont. from prev. page)

STRATAGENE EAGLE EYE II 04/07/98 18:54:46

IMAGE SIZE (640 x 480 x 8).
INT PERIOD = 0.19 SEC.
ACQUIRED TUE APR 07 18:53:02 1998.

loaded 10 pl / lane, the same gel

all the bands are too faint for box traces

10x Terg poly buffer 150

dNTPs (25mM) 12

EC1605 1

CC125 1

Taq 2000 10

dH₂O 1326

added 47.5 µl ea reaction

95°C, 3' / 55°C, 3' / 72°C, 2' // 1 cycle

95°C, 1' / 55°C, 1' / 72°C, 2' // 24 cycles

72°C, 10' → 6°C, ∅


Confirmation digests of RIL fragment used for LIC cloning of pACYC-RIL and mini-prep of selected ESP4Flp, ESP4Kw80 and ESP4Kw70 clonesextracted fragment sizes:

STRATAGENE EAGLE EYE II 04/07/98 18:58:02

IMAGE SIZE (640 x 480 x 8).
INT PERIOD = 0.33 SEC.
ACQUIRED TUE APR 07 18:48:13 1998.

ESP4Flp ESP4Kw80 ESP4Kw70

NotI EcoRV #1

332
ESP4Kw70: EcoRV: 7727, 2569, 1660HindIII: 5841, 2398, 2300, 1022
(not separated)

ESP4Kw80 EcoRV: 7727, 2549, 1660

ESP4Flp NotI: 9111, 1857, 1898

(maps of ESP2Kw70 and ESP2Kw80 next page)

all digest result in the expected pattern

test of RIL fragment

undigested: 610 bp

NotI: 220 bp, 390 bp (w/ salt retards smilling.)

XbaI: 180, 430 bp (w/ salt retards)

RIL fragment seems to be correct

under why I can't clone to toxic?

STRATAGENE IN-HOUSE OLIGONUCLEOTIDE REQUEST FORM

REQUESTED BY: Carten Carten DATE: 4/22/91
 ACCOUNT #S: (BUSINESS) 12-33050 (PERSONAL) 12-12

SEQUENCE(S) 5' → 3':

1. Argu fund / SpeI
 5'-GAGIAGTIAAGTAAATCAGIACGICGGTTCG-
 TCGIAG- / - - - - -
 27 bp
2. Low W rev / SpeI
 5'-GAGICTGAGAGIGGGIATGCGATCAIACG-
 CTATGCA- / - - - - -
 30 bp
3. HeY rev / SpeI
 5'-GAGTCTGAGTICATGAGTGTITATTC-
 CGTGG- / - - - - -
 27 bp
4. HeY fund / SpeI
 5'-GAGIAGTIAAGTICGTTGAIATIGGCIGTT-
 AGTICA- / - - - - -
 27 bp

5' MODIFICATIONS: 1. _____ 2. _____ 3. _____ 4. _____

PURIFICATION METHOD (circle one for each):

1. Etoh ppt
PAGE

2. Etoh ppt
PAGE

3. Etoh ppt
PAGE

4. Etoh ppt
PAGE

COMMENTS/SPECIAL REQUESTS:

(oligo department use only)

DATE/TIME REC'D: _____ DATE NEEDED: _____
 COMPLETED BY: _____ DATE COMPLETED: _____

ordered new primer
for construction of RL
array. LIC cloning
appeared to be very
ineffective in previous
attempts.

→ try to clone between
the SpeI and KpnI sites
of pACYC-tet

also ordered forward
and reverse primer for
RI and IL array
to be cloned the same way

Test of production of Ku80 and Atp as GST fusion proteins in yeast

goal: produce the said proteins as GST fusion proteins. Previously,
failed as CBP fusion proteins either in E. coli or yeast as
a host

picked colonies from SP-Q01 yeast cells transformed
with ESP4 Ku80 and ESP4 Atp (pages 68 and 70) on
EtoH plates → incubated at 30°C

C.C.D.

(Continued)

↓

Split each culture into 2x 5ml cultures

↓

pelleted cells from each 5ml culture (10' at 2500 rpm, tabletop)

↓

resuspended 2x each 50ml sterile ddH₂O

↓

resuspended each pellet in 10ml unsequenced EAT

Added to one culture each group 10 µl thiocyanate (5 mM)

Becomes "repressed/uninduced" sample


↓

incubated at 30°C in shaker at 250 rpm (6h incubation)

PCR amplification of RIL, IL and RI fragments using primers with
 SpeI/XhoI extensions

GENE EAGLE EYE II 04/30/98 11:39:11

148 1998.



	1	2	3
DxTaq + Pcc.	10	10	10
dUTPs	.8	.8	.8
IL-RIL (page 47)	1	1	1
CC 972B (tag/ufwd)	.5	.5	-
CC 973 B (leuW rev.)	.5	-	.5
CC 974 B (ileY rev.)	-	.5	-
CC 975 B (ileY fwd)	-	-	.5
Taq + Pcc.	1	1	1
ddH ₂ O	86.2	86.2	86.2

95°C, 3' / 55°C, 3' / 72°C, 2' // 1 cycle

95°C, 1' / 55°C, 1' / 72°C, 2' // 19 cycles

72°C, 10' → 6°C hold

isolated PCR products using Stratagene PCR purification kit.

eluted DNA into 50 µl 1x Universal buffer

- added 2 µl SpeI and 2 µl XhoI to each purified reaction

- incubated 29 at 37°C

- loaded on 1% 1.5% agarose gel

(Continued)

C.O.S.D.

(Continued)

↓
1 gel transferred to nitrocellulose membrane
by Southern transfer

↓
kept in TBS, 5% dry milk at 4°C in cold room

other gel stained with Coomassie

construction of pACYC-RIC, pACYC-RI and pACYC-IL

a) beaker digest

10 μ l pACYC-RIC miniprep

5 μ l 10 \times U.B.

1 μ l SpeI

1 μ l XhoI

83 μ l ddH₂O

incubated 1 h at 37°C

removed enzymes with ^{beaker to ready} PCR purification buffer (Stratagene
PCR purification kit). Purified through spin cups and
eluted into 50 μ l 20 mM Tris-HCl pH 8.0

b) ligation

2 μ l beaker (pACYC-LIC)

5 μ l RIC, RI or IL (page 84/85)

4 μ l 5 \times ligation buffer

1 μ l T4 DNA ligase

8 μ l ddH₂O

incubated for 4 h at 16°C (rotorcycle cold block)

↓
used 2 μ l to transform XL1 blue


C. C.

PCR screen of pACYC-RIL and pACYC-IL (see prev. page)

	A	B	
10x Targ Poly buffer	100	100	95°C, 3' / 75°C, 3' / 72°C, 2' // 1 cycle
dUTP (25 mM sol.)	8	8	95°C, 1' / 55°C, 1' / 72°C, 2' // 29 cycles
CC 972 B (argU fwd)	1	1	72°C, 10' → 6°C, 5 min
CC 973 B (LeuW rev.)	1	-	and cut up
CC 974 B (LeuY rev.)	-	1	
Targ 2000	10	10	loaded 10 µl reaction on 1.2%
ddH ₂ O	880	880	agarose gel

STRATAGENE EAGLE EYE II 05/04/98 16:51:48

IMAGE SIZE (640 x 480 x 8).
INT PERIOD = 0.09 SEC.
ACQUIRED MON MAY 04 16:49:57 1998.



Expected sizes:

RIL 610 bp

IL 430 bp

pACYC-RIL: clone #7


pACYC-RIL: 23, 25, 28, 29, 32, 33

Coomassie gel of induction of P_{leu} and K_u80 in yeast (contin. of page 95)

STRATAGENE EAGLE EYE II 05/04/98 15:38:24

FILE C:\EAGLE\YE\FLP80SP2.TIF
IMAGE SIZE (640 x 480 x 8).
REAL-TIME ACQUIRE.

ACQUIRED MON MAY 04 15:36:15 1998



lanes slightly overloaded
only the GST control (ES12)
has an induced band visible
(arrow)

there may be a very weak
induced band at ~110 kDa in
the K_u80 sample, which
would be roughly the
expected size

C. Co

(Continued)

digest vector preps for 2h at 37°C


↓
 isolated DNA using Strataprep PCR purification kit
 (150 µl DNA concn. soln. → on spin cups)
 eluted DNA into 50 µl TE

ligation:
 2 µl vector preps (C2 or C6)
 5 µl inserts (L1A or L1B)
 4 µl 5x ligase buffer
 1 µl T4 DNA ligase
 8 µl dd H₂O
 ligated at 4°C (overnight) o/n

STRATAGENE EAGLE EYE II 05/06/98 17:02:28

 IZE (640 × 480 × 8).
 IOD = 0.33 SEC.

0 WED MAY 06 17:00:54 1998.



manipulations of pACYC-R1 and
 pACYC-RIL colonies (continuation
 page 87)

expected fragment sizes

2192 bp, 1230 (R1) or ⁴⁰⁰~~1230~~ (RIL)

both pACYC-R1 clones
 are confirmed by restriction
 digests.

None of the RIL clones is
 positive although they were
 PCR positive (comp. page 87)


C. CO

re-screened pATC-RIL colonies by PCR

10x Terg - poly 12.5/1hr	200	94°C, 3' / 55°C, 3' / 72°C, 2' // 1 cycle
oligos	16	94°C, 1' / 55°C, 1' / 72°C, 2' // 29 cycles
CC 972 B	2	72°C, 10' → 6°C, hold
CC 973 B	2	used last top
Tag 2000	10	
dH ₂ O	1770	loaded 5µl on 1% agarose gel

STRATAGENE EAGLE EYE II 05/06/98 16:02:19

IMAGE SIZE (640 x 480 x 8).
INT PERIOD = 0.16 SEC.
ACQUIRED WED MAY 06 16:00:05 1998.




None of the previous
pATC-RIL colonies were
positive.

3 of the new colonies
(49, 50 and 55) are
positive

STRATAGENE EAGLE EYE II 05/07/98 14:57:3

IMAGE SIZE (640 x 480 x 8).
INT PERIOD = 0.49 SEC.
ACQUIRED THU MAY 07 14:55:05 1998.



Miniprep of pATC-RIL clones
(above)

expected

Not: 2192 bp, 1400 bp

YbATc

but all clones are positive by
NotI cuts but YbATc site
seems to be missing
(either from vector or from insert)

C. C. S.

total Fbp expression in pACTC-RIL/BL21DE3 cells

picked colonies from transformations on page 93 (one each expression clone in each RIL strain) and inoculated o/n cultures (1ml Amp/can LB)

induction of Fbp expression in BL21DE3/RIL strains

inoculated 1ml LB (no antibiotic) with 30 μ l of o/n culture

↓
grew for 2.5h at 37°C in shaker

↓
removed 100 μ l culture, added loading buffer and boiled for 2' → uninduced samples

↓
added 2 μ l IPTG (500mM → 1mM final) to pCellFbp samples

↓
incubated for another 2h at 37°C in shaker

↓
removed 100 μ l from culture, added 100 μ l SDS loading buffer and boiled for 2' → induced samples

↓
let uninduced E. coli grow for another 2.5h (induced)

↓
removed 100 μ l culture, added 100 μ l SDS loading buffer and boiled for 2' → induced samples

Next attempt at LIC cloning Fbp into pACTC-LIC

4 μ l LIC-tailed Fbp fragment (page 81)

1 μ l LIC-tailed pACTC-LIC vector (from Quimble)

5 μ l dH₂O

incubated 1h at RT, extracted 3x with Stratamclean

→ transformed XL1-blue with 1 μ l of LIC reaction → on LB Cam

CCB

Detection of F10 expression in pRBL21 DE3 RIL strains

loaded 2 lanes 4-20% gels each with pColF10₁ induced/uninduced
in RIL 49, 50 and 55 (only induced) and E. coli F10 induced/uninduced
in RIL 49 and 55

run gels 1.5h at 150V

one gel stained with coomassie

other gel transferred to nitrocellulose membrane by
Sambrook transfer (Nover apparatus) for 2h at 30mA

Western blotting of F10 clones with new anti CBP antibody from
antibody production group.

blocked membrane with 5% dry milk in TBS for 1h

incubated membrane with 10ml TBS, 5% dry milk
rabbit anti CBP polyclonal 1:1000 dilution for 1h

washed 3x 5' with TBST

incubated 1h with 10ml HRP-conjugated anti rabbit AB
in TBS/5% dry milk

washed 3x 5' with TBST
1x 5' with TBS

developed with Pierce Supersignal system

Signal extremely strong. Had to add piece of paper
between the membrane and the film to achieve reasonable
exposure times (chromoluminogram next page).
took image with Eagle eye (signal strong enough for that)

C. O.